PTB PCA

November 23, 2021

0.0.1 Read MIT format .dat ecg data files and .hea headers

```
[8]: #Download first at: https://physionet.org/content/ptbdb/1.0.0/
     #BASE DIR = '/media/julian/Volume/data/ECG/mit-bih-arrhythmia-database-1.0.0/'
     →#Arrhythmia
     BASE_DIR = '/media/julian/Volume/data/ECG/ptb-diagnostic-ecg-database-1.0.0/'
     def get_file_list(BASE_DIR, relative=True):
         record_files = []
         \#file\_endings = ['.dat', '.hea', '.xyz']
         with open(os.path.join(BASE_DIR, 'RECORDS')) as recs:
             record_files = recs.read().splitlines()
         if not relative:
             record_files = [os.path.join(BASE_DIR, f) for f in record_files]
         return record files
     record_files = get_file_list(BASE_DIR)
     print(len(record_files))
     print(record_files[:10])
    549
```

```
['patient001/s0010_re', 'patient001/s0014lre', 'patient001/s0016lre', 'patient002/s0015lre', 'patient003/s0017lre', 'patient004/s0020are', 'patient004/s0020bre', 'patient005/s0021are', 'patient005/s0025lre']
```


0.0.2 Extract signal from *.dat files & Read annotations & Read comments


```
[9]: def plot_cov(co, title=None):
          if title:
              plt.title(title)
          plt.imshow(co)
          plt.colorbar()
          plt.show()
[10]: def read_comment_map_PTB(record_path):
          #print(record_path)
          record = wfdb.rdrecord(record_path)
          comment_map = {}
          for c in record.comments:
              e = c.split(':')
              comment_map[e[0]] = e[1].strip()
          return comment_map
[11]: def filter_comment(comment, key):
          c = comment
          if key == 'Reason for admission':
              if 'Cardiomyopathy' in c or 'Heart failure' in c:
                  return 'Cardiomyopathy'
              elif 'n/a' in k or 'Palpitation' in k:
                  return 'Miscellaneous'
              elif 'angina' in k:
                  new_comments['Angina'] = comments[k]
              else:
                  new_comments[k] = comments[k]
[12]: def read_comment(record_path):
          record = wfdb.rdrecord(record path)
          return record.comments
[13]: def read_header(record_path):
          record = wfdb.rdheader(record_path, rd_segments=True)
          return record.comments
[14]: def read_signal(record_path, physical=True):
          #print(record_path)
          record = wfdb.rdrecord(record_path, physical=physical)
          #print_object_attributes(record)
          if physical:
              data = record.p_signal
              data = record.d_signal
          return data
```


Save all signals and attributes in file_data (also note how many had functioning annotations)


```
[16]: file_data = []
success = 0
for f in record_files[:]:
    p = os.path.join(BASE_DIR, f)
    d = read_signal(p, physical=False)
    file_data.append((f, d))
```


```
[18]: import numpy as np
  from numpy import linalg as LA
  from sklearn.preprocessing import normalize
  for f in file_data:
        data = normalize(f[1], norm='12').T
        co = np.cov(data)
        w,v = LA.eig(co)
        w = np.expand_dims(w, 0)
        plot_cov(co, f[0])
        plot_cov(w, f[0])
```




















[]: